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en	INSTRUCTIONS FOR USE – INFUSION SETS
<b>Intended Use</b>	Standalone devices for the connection to, and the extension of, intravenous infusion systems, incorporating valves and other features for use during the infusion of fluids into the human body.
<b>Indication(s)</b>	<b>A-Sets:</b> Patients that require pumped drug infusions (syringe driver or volumetric pump), such as in Total Intravenous Anaesthesia (TIVA). <b>Safe-Sets:</b> Patients that require intermittent drug administration, drug and fluid administration remote from the patient's cannula
<b>Contraindications</b>	Do not use on a patient with known hypersensitivity to silicone or plastic materials. Please contact the manufacturer for advice if in doubt.
<b>Sterility</b>	Sterilised by Ethylene Oxide, do not use if pouch is damaged
<b>Intended Patient Group(s)</b>	No restrictions associated with any patient characteristics such as age, gender, race or pre-existing medical conditions such as pregnancy.
<b>Intended User(s) &amp; Facilities</b>	The devices are intended to be used by health professionals only, within a controlled clinical environment according to national regulations.
<b>Clinical Benefits &amp; Performance Characteristics</b>	The clinical benefits are typically derived from the infused substances rather than the infusion sets themselves, however the sets have been designed with appropriate characteristics to ensure infusion is performed reliably and safely, i.e.: Luer lock connections, anti-syphon valve on drug line(s), anti-reflux valve on fluid line(s) and minimal dead space where drug & fluid lines join.
<b>Warnings / Precautions for Storage, Handling, Preparation &amp; Use</b>	Ensure the devices are stored in a dry environment away from direct sunlight. Make sure infusion lines are set up to avoid the possibility of kinks forming during use. Use clamps to occlude pressure line(s) when connecting to the cannula and when it is not needed for infusion purposes. Only use with compatible Luer lock devices (e.g., syringe, cannula...etc.). Swab needle-free valves before use. <b>Ensure all lines are primed before use.</b> Do not reuse the device as this may cause cross infections and/or drug reactions. Dispose via clinical waste streams – used sets would pollute recycling streams.
<b>Residual Risks &amp; Undesirable Side-Effects</b>	If not used appropriately according to the instructions above, there is a risk of: <ul style="list-style-type: none"> <li>• Occlusion or backflow causing under-infusion or over-infusion</li> <li>• Syphoning from pump or syringe causing accidental overdose</li> <li>• Patient drug reactions from premature or inappropriate mixing</li> <li>• Patient allergic reaction if hypersensitive to materials</li> <li>• Embolism from admission of air into administration lines</li> </ul>
<b>In the event of an incident or defective device</b>	If any serious incident has occurred in relation to the device, the user and/or patient should be report it to the manufacturer via <a href="mailto:info@amedicalsystems.co.uk">info@amedicalsystems.co.uk</a> , and the competent authority of the Member State in which the user and/or patient is established (refer to <a href="https://ec.europa.eu/health/md_sector/contact_en">https://ec.europa.eu/health/md_sector/contact_en</a> )

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